



AED.5

Comparison of our AEDs



Comparison	ZOLL AED3	STRYKER SAMARITAN	PHILIPS Onsite	STRYKER CR2
Available in fully automatic mode	✓	✓	✗	✓
Available in semi-automatic mode	✓	✓	✓	✓
Electrode life:	5 years	4 years	2 years	4 years
Universal Adult / Child electrodes	✓	✗	✗	✓
Feedback on the depth and speed of compressions	✓	✓	✗	✓
Guarantee	6 years + 2 years with registration	8 years	8 years	8 years
Technical details				
AED management program via NON WIFI - Alerts and notifications on computers or mobiles	✓	✗	✗	✓
Alerts and notifications on NOT computers or mobiles when the AED is in motion and / or in use	✗	✗	✗	✓
Help / Guide / Guidelines for CPR	✓	✓	✗	✓
Guide / Metronome to guide the compression speed	✓	✓	✓	✓
Heart rate analysis at the same time as CPR	✓	✗	✗	✓
Background noise analysis and volume correction	✗	✗	✗	✓
Detects and abstracts surrounding vibrations	✗	✓	✗	✓

Comparison	ZOLL AED3	STRYKER SAMARITAN	PHILIPS Onsite	STRYKER CR2
Calculation of body impedance	✓	✓	✓	✓
Shock power:	120 to 200 joules	150 to 200 joules	150 joules	200 to 360 joules
Shock charging time:	8 seconds (semi-automatic model), 13 seconds (automatic model)	8 seconds	8 seconds	0 seconds, Charged shock during CPR
Sending live event data to doctors and paramedics	X	X	X	✓
Electrodes without packaging	✓	✓	✓	✓
Button to switch to child mode	✓	X	X	✓
Intervention kit included (CPR mask, razor, wipe, gloves, scissors cut everything)	✓	✓	✓	✓
Battery life:	5 years	4 years	5 years	4 years
Language change in 1 English / French button	X	X	X	✓
Recording of event data	✓	✓	✓	✓
Transfer to computer and access to event data	✓	X	X	✓
Size weight :	12.7 x 23.6 x 24.7 cm, 5,5 lbs	20 x 18,4 x 4,8 cm, 2,4 lbs	7 x 19 x 21 cm, 3,3 lbs	9,7 x 22,6 x 27,4 cm, 4,5 lbs
Protection index against dust and water:	IP55	IP56	IP2X	IP55